### **DATA SHEET**

# VETRAP 3M DOSAGE AND METHOD OF USE

Under no circumstances should the band be applied with a higher tension than 50% of the maximum tension.

### SPECIAL WARNINGS

Vetrap 5 cm: blue, red, greenVetrap 7,5 cm: blue, red, green

• Vetrap 10 cm: blue, fuchsia, yellow, black, green, red



#### **PACKAGING**

- Vetrap 5 cm:box of 200 pieces (can be sold individually)
- **Vetrap 7,5 cm:** box of 150 pieces (can be sold individually)
- Vetrap 10 cm: box of 100 pieces (can be sold individually)

## **DESCRIPTION**

Elastic bandage for restraint and maintenance. Lightweight, strong and porous, suitable for wound protection, holding dressings in place and tail dressing.

- Adheres to itself without the need for additional fixings
- maintains mobility
- easy to remove, does not stick to the hair
- allows the bandaged area to breathe
- water- and sweat-proof
- if sufficient pressure is applied, it can be used to stop bleeding while remaining flexible enough not to impair blood circulation

