

## DATA SHEET

# FLEX-1-ONE

## DOSAGE AND METHOD OF USE

Under no circumstances should the band be applied with a tension greater than 50% of the maximum tension

### SPECIAL WARNINGS

- **Flex-1-One 5 cm:** blue, yellow, red, green
- **Flex-1-One 7,5 cm:** white, blue, fuchsia, yellow, red, green
- **Flex-1-One 10 cm:** white, blue, fuchsia, yellow, black, red, green



### PACKAGING

- **Flex-1-One 5 cm:** box containing 24 bands
- **Flex-1-One 7,5 cm:** box containing 12 bands
- **Flex-1-One 10 cm:** box containing 12 bands

### DESCRIPTION

Cohesive bandage useful for the maintenance of bandages, post-operative dressings and tail dressings.

Its characteristics are:

- high resistance to rubbing
- is easily torn by hand without the need for scissors,
- controlled compression (not constrictive)
- adheres to itself without being sticky on the hair
- water- and sweat-proof
- provides support for tendons and ligaments, maintaining position even in particularly difficult conditions

